



## EZ-SCAN<sup>tm</sup>

### A new, compact, bench-top TLC scanner / detector system

The EZ-SCAN system occupies only a small footprint on scarce lab-bench real-estate. Using EZ-SCAN with the Omni-Rad<sup>®</sup> detector module you can scan a small sample drawn directly from a production run, without requiring any dilution.

The annotated thumb-nail image shows a 'mock' scan on a TLC plate on which small drops of <sup>18</sup>F-labeled fluoride were placed by hand. The plate was later cut up and the pieces assayed.

